

900RD110 Humidified Resuscitation System

User Instructions


①


② ③ ④

⑤

⑥

900RD110 Humidified Resuscitation System

User Instructions

Set Up

Consult the accompanying documents for your Neopuff™ Infant T-Piece Resuscitation Device, MR850/MR810 Humidifier and MR290/MR225 Chambers for additional user instructions and warnings.

① Fit a Chamber to the Humidifier

Slide a MR290 chamber onto the humidifier base. Remove the blue caps. Hang the water bag, unwind the waterfeed set and spike water bag. The bag should be at least 50cm above the chamber. Use USP sterile water for inhalation or equivalent.

Or

Slide a MR225 chamber onto the humidifier base. Remove the blue caps. Fill with 30mls of sterile water or equivalent. This will allow for 60 minutes of running time (@ 8 L/min flow rate) and then refill as necessary during resuscitation.

② Connect Circuit

- a) Connect the dry-line between the chamber inlet port and the Neopuff™ gas outlet.
- b) Connect the blue inspiratory limb to the remaining chamber port.
- c) Ensure the temperature probes and heater wire adaptors are fully inserted. If temperature probes are not used (i.e. using the MR810), ensure caps are placed over the temperature probe ports.

Connect to Flow Source

- ③ a) Connect an oxygen or blended oxygen/air supply to gas inlet port using gas supply line.
- b) Adjust gas supply to desired flowrate between 8 and 15 L/min. Recommended operating gas flow rate: 8 L/min.

Set Humidifier

- ④ Turn on the humidifier base.
- MR850 – Check the humidifier is set to invasive mode. if not hold the mode button until the invasive mode light is lit.
- MR810 – Set the humidifier to high mode. Ensure the green heater wire light is lit. Caution is required when delivering gases that are not body temperature saturated, to patients with bypassed airways.

Set PIP

- ⑤ Ensure the test lung is inserted into the patient T-Piece. Occlude PEEP (Positive End Expiratory Pressure) cap for adjustment and setting of PIP (Peak Inspiratory Pressure).

⑥ Set PEEP

Adjustment of PEEP.

To Resuscitate

See RD900 Neopuff™ Infant Resuscitator Operating Manual.

- Connect the test lung to the T-Piece circuit and test resuscitator function before connecting to a patient.

Intended Use

The Fisher & Paykel Humidified Resuscitation System is designed to be used for Infant Resuscitation to deliver appropriate PEEP and PIP Pressures. This circuit (900RD110) is designed for use with the F&P Neopuff™ Infant Resuscitator and F&P humidifiers only.

Product Specifications

Input Flow Range	8 – 15 L/min	
Circuit Length	Inspiratory	1.3m
	Dry Line	0.6m
Positive End-Expiratory (PEEP)	@8 L/min 1 to 10 cmH ₂ O [mbar] @10 L/min 2 to 13 cmH ₂ O [mbar] @15 L/min 4 to 14 cmH ₂ O [mbar]	
Device Connection	10 mm Female	
T-Piece Connection	15 mm Female (Conical Connector)	
Disposal	Discard according to hospital protocol	
Environmental limits	Operating: 18 °C to 26 °C (+64 to +79 °F), 30 to 75% humidity Storage: -10 °C to 50 °C (+14 to +122 °F), Up to 95% humidity	
Patient Mass Range	<10 kg (22 lb)	
Humidity Performance	MR850	>33 mg/L
	MR810	>10 mg/L

Warnings

- This product is intended to be used for a maximum of 7 days.
- This circuit is for use during resuscitation ONLY.
- DO NOT cover the circuit with materials such as blankets, towels or bed linen.
- Avoid prolonged contact with patient skin.
- DO NOT stretch or milk tubing.

The Humidified Resuscitation system is designed for SINGLE PATIENT USE ONLY. Reuse may result in transmission of infectious substances, interruption to treatment, serious harm or death.

DO NOT soak, wash, sterilize or re-use this product.

DO NOT occlude the resuscitation tubing.

Ensure no smoking, naked flames or sources of ignition are present when using this circuit.

The Humidified Resuscitation system is only to be used with a flow rate from 8-15 L/min. DO NOT attempt to use product outside this flow range.

The T-Piece circuit Humidified system and Neopuff™ are intended for use by adequately trained medical practitioners.

Ensure pressures are monitored throughout resuscitation.

Rx only.

DO NOT operate the humidifier or the circuit with no gas flow.

Ensure the humidifier is mounted below the level of the patient.

The Humidified Resuscitation System is to be used with Medical Grade oxygen/air.

DO NOT fill the chamber above the maximum fill line.

Latex Free

Do Not Use If Package Is Damaged

California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

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REF 185043544 REV J 2019-04 © 2019 Fisher & Paykel Healthcare Limited