
F&P FlexiTrunk™ Nasal CPAP Interface User Instructions

Manufacturer ■ Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo - SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132 **Denmark** Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France** **REP** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D 73614 Schorndorf, Germany Tel: +49 7181 98599 0 Fax: +49 7181 98599 66 **Hong Kong** Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115 **Korea** Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: +52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 00800 49 119 77 Fax: 00800 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50 **Spain** Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Ostim Mahallesi 1249, Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001

F&P and FlexiTrunk are trademarks of Fisher & Paykel Healthcare Limited.
For patent information, see www.fphcare.com/ip

REF 185042952 REV W 2021-03 © 2021 Fisher & Paykel Healthcare Limited

Rx only 0123

Fisher & Paykel
HEALTHCARE

FlexiTrunk™ Nasal CPAP Interface

User Instructions

en

Intended Use

The FlexiTrunk nasal CPAP Interface is designed to deliver non-invasive positive pressure therapy to a spontaneously breathing patient weighing up to 10 kg in a hospital or clinical setting where the patient is adequately monitored by trained medical staff.

Product Specifications

BC190 BC191 BC192	The FlexiTrunk nasal CPAP Interface is designed to connect to the Fisher & Paykel Healthcare Bubble CPAP System. The adaptors allow connection to other manufacturers' CPAP circuits with 7.8 or 10 mm connectors.
Maximum input flow	15 L/min
Maximum allowable operating pressure	15 cmH ₂ O

Please refer to last page for dimensions and components.

Clinical Benefit

It is established in the clinical evidence that CPAP is a form of non-invasive respiratory support which can reduce the frequency and duration of invasive mechanical ventilation and the associated complications of intubation and mechanical lung injury.

Contraindications

The following are contraindications related to nasal CPAP and the use of nasal CPAP interfaces:

- Non-spontaneous breathing
- Congenital abnormalities or malformations where bi-nasal prongs or nasal mask are contraindicated (e.g. choanal atresia).
- Congenital abnormalities or malformations where positive pressure therapies are contraindicated (e.g. diaphragmatic hernia and tracheo-oesophageal fistula).
- Nasal trauma/severe deformity that might be exacerbated by use of nasal prongs or nasal mask.

Side effects

The following are the known side effects related to nasal CPAP and the use of nasal CPAP interfaces:

- Gastric distension and feeding intolerance of minor severity are undesirable side effects of using CPAP therapy.
- Nasal trauma and skin injury of minor to moderate severity are undesirable side effects of using CPAP interfaces.

WARNINGS

- The components of the FlexiTrunk Interface are designed for single patient use only. Reuse may result in transmission of infectious substances, interruption to treatment, serious injury or death.
- Handle with care. Use caution when positioning or disconnecting the Interface. Avoid excessive pull forces, sharp objects and tubing holders. Damage to the tubing may cause loss of pressure and require immediate replacement.
- This product is intended to be used for a maximum of 7 days.
- DO NOT use if product or packaging is damaged.
- US Federal law restricts this device to sale by, or on the order of, a physician.
- DO NOT modify this product.
- Use patient oxygen monitoring.
- Connection to other manufacturers' humidification systems increases risk of condensate and/or high gas temperature delivery to the patient.
- Connect only pressure monitoring devices to the pressure monitoring port on the nasal tubing.
- FlexiTrunk Interface components are supplied in a clean state. Employ suitable methods to prevent bio-contamination during use and upon disposal.
- DO NOT soak, wash, sterilize, or reuse this product.
- DO NOT overtighten the bonnet straps.
- DO NOT cover the patient's face with the bonnet as suffocation may occur.
- USE ONLY the recommended FlexiTrunk Interface components, impaired positioning of the nasal prongs and mask can occur from non-approved substitutes.
- Septal Necrosis can occur if correct sizing or placement of prongs and mask is not observed.
- Always use pressure monitoring to verify that the patient is receiving the prescribed CPAP level.
- DO NOT use medications containing Tyloxapal (such as Tacholiquin) as this may damage the tubing and lead to loss of CPAP pressure.
- The use of this device is not without risk. Even if used as intended following all instructions and warnings provided, the following risks remain: airway injury, hypoxic injury, infection, mechanical lung injury, skin burns, nasal trauma and skin injury. These risks may result in serious injury or death.
- California residents please be advised of the following, pursuant to Proposition 65: This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65
- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.

Symbol Definitions

	Consult Instructions for use	CE 0123	CE Mark		Not made with phthalates (DEHP, DBP, BBP)
	Latex free		Do not reuse	REF	Catalogue number
EC REP	Authorized representative in the European Community	Rx only	Prescription only		Manufacturer
	Date of manufacture		Use by date		Recyclable PET symbol
	Storage temperature limit	MD	Medical device	LOT	Batch code
CE	CE Mark		Caution		Handle with care
S	Size Small	M	Size Medium	L	Size Large
XL	Size Extra Large				

Set-up

- Refer to the F&P Bubble CPAP System instructions prior to set up of the FlexiTrunk Interface.
- If using the adaptors, refer to the set-up instructions provided by the circuit manufacturer prior to set up of the FlexiTrunk Interface.
- Refer to illustrations 1-15 for correct set up.
- For correct product sizing, refer to the tables on the last page.
- A maximum pressure relief device MUST be used. Pressure relief settings shall be suitable for infant and neonatal CPAP.

Checks Before Operation

- When pressure monitor port on nasal tubing is not in use, the cap must be closed.
- Check that all circuit connections are tight before use and after any adjustment.

Set-up Instructions


Refer to illustrations on page 1

- ① Measure the baby's head circumference in centimeters.
- ② Choose the correct sized bonnet. The bonnet should be snug fitting. Stretch the bonnet with your hands initially for ease of placement.
- ③ Slip the bonnet onto the baby's head, completely covering the ears and with the back edge of the bonnet at the base of the neck. The front edge of the bonnet should be just above or on the eyebrows.
- ④ Use the sizing guide to choose suitable prongs or mask. Prongs should fill the nares completely without stretching the skin. Use the biggest possible size. Mask should not touch the edge of the nose, septum or the eyes.
- ⑤ Choose the correct length of the nasal tubing. Use the shortest length possible. The clear tubing should not extend over the patient's forehead. The following may be used as a guide:
50 mm ≤ 1.5 kg, 70 mm ≤ 2.5 kg, 100 mm > 2.5 kg
- ⑥ Connect prongs or mask to nasal tubing ensuring that it is inserted fully. If using prongs: squeeze sides of prongs firmly to expose grooves. Start from one end and insert prong grooves into nasal tubing. Push end in firmly.
- ⑦ Remove tear-off strips as required to adjust the nasal tubing angle to optimise the seal. DO NOT discard tear-off strips as they may be needed when alternating between nasal prongs and masks.
- ⑧ Attach the circuit to the nasal tubing in either orientation. Refer to instructions provided by circuit manufacturer for more detail. Note: Disconnect FlexiTrunk gently by twisting the connectors. DO NOT pull forcefully, as this may damage the tubing. Handle with care.
- ⑨ Set gas flow to prescribed level. Turn the humidifier on. Refer to manufacturer's instructions. **Place hand close to nasal prongs or mask to ensure there is gas flow present.**
- ⑩ If using prongs: Clear nasal secretions before inserting the nasal prongs. **Ensure the nasal prongs are positioned at least 2 mm (0.08 inches) from the septum to avoid pressure necrosis. Adjust as necessary. Hourly checks of the septum integrity are recommended.**
- ⑪ If using mask: Connect to patient by placing mask around nose. The mask should sit comfortably around the patient's nose. It must not occlude the nostrils or touch the septum and should not be over the lip or over the eyes.
- ⑫ Set height of the foam block by adding or removing tear-off strips. Setups with masks may require adding strips for proper positioning. Place the strap over the foam block and tighten to hold the tubing firmly in place.
- ⑬ Hook the clips from the side straps of the bonnet to the glider. Pull both straps at the same time for central positioning. Affix the round velcro tabs to the blue strap of the bonnet. Use the least tension possible to maintain prescribed CPAP level and stability of the Interface.


- ⑭ DO NOT use angel frames, Christmas trees or similar devices as they can puncture the FlexiTrunk tubing and cause loss of pressure.
- ⑮ Ensure the weight of the circuit is gently supported to reduce tension on the nasal tubing. A pillow or folded towel is recommended. The FlexiTrunk tubing must have free movement to allow the baby's head to move without restriction. The circuit should be kept below the level of the nasal tubing to minimize condensate at the baby.


During Use:


- Check patient frequently for signs of redness, pressure sores or irritation which may result from extended prong or mask use. Hourly checks are recommended.
- Alternating between mask and prongs can reduce the risk of irritation. Alternating every 4 to 6 hours is recommended.
- Positioning of the mask should be adjusted for best fit and minimal leak. If excessive leak is still present, flow may be increased to compensate. Maximum allowable flow rate is 15 L/min.
- Check condensate regularly and drain as required.
- Replace the bonnet if it becomes loose during prolonged usage.

FlexiTrunk Nasal Prongs		
		
Product Code	A	B
BC3020	3.0 mm	2.0 mm
BC3520	3.5 mm	2.0 mm
BC4030	4.0 mm	3.0 mm
BC4540	4.5 mm	4.0 mm
BC5040	5.0 mm	4.0 mm
BC5050	5.0 mm	5.0 mm
BC5550	5.5 mm	5.0 mm
BC5560	5.5 mm	6.0 mm
BC6060	6.0 mm	6.0 mm
BC6070	6.0 mm	7.0 mm
BC6570	6.5 mm	7.0 mm

FlexiTrunk Nasal Tubing	
	
BC19x series	
Product Code	Size
BC190	50 mm
BC191	70 mm
BC192	100 mm

Chinstrap	
	
Product Code	Size
BC351	20–26 cm
BC353	26–32 cm
BC355	32–38 cm
BC357	38–44 cm

FlexiTrunk Nasal Mask	
	
Product Code	Size
BC800	S
BC801	M
BC802	L
BC803	XL

FlexiTrunk Bonnet	
	
Product Code	Size
BC300	17–22 cm
BC303	22–25 cm
BC306	25–29 cm
BC309	29–36 cm

FlexiTrunk Headgear	
	
Product Code	Size
BC325	29–36 cm
BC328	35–40 cm
BC331	40–45 cm

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