





RD1300-10 T-Piece Circuit (with duckbill port)

User Instructions

RD1300-10 T-Piece Circuit (with duckbill port)

User Instructions

Set-up

⚠️ Consult the accompanying documents for your T-Piece Resuscitation Device for additional user instructions and warnings.

⚠️ An alternative means of resuscitation must be available

- ① Connect the RD1300-10 to the T-Piece Resuscitation Device.
 - Ensure device is connected to recommended flow source
 - Adjust flow to desired flow between 5 and 15 L/min
- ② Occlude PEEP (Positive End Expiratory Pressure) cap for adjustment and setting of PIP (Peak Inspiratory Pressure)
- ③ Adjustment of PEEP

To Resuscitate

⚠️ Test resuscitator function with blue cap in place on the T-Piece circuit or connect test lung to T-Piece circuit before connection to the patient.

⚠️ Ensure you have removed the blue cap or test lung before connecting to the mask or endotracheal tube.

⚠️ Ensure pressures are monitored throughout resuscitation.


- ④ Inspiration
- ⑤ Expiration

Duckbill Port

The duckbill port accepts catheters from 5 Fr to 8 Fr (1.7-2.7 mm OD, 0.06"-0.11" OD)

- ⑥ Delivering surfactant through port
- ⑦ Suctioning through port

Rx only

Intended Use

The T-Piece Circuit (RD1300-10) is designed to be used for infant resuscitation to deliver appropriate PEEP and PIP pressures when used with a flowrate from 5-15 L/min. The T-Piece Circuit (RD1300-10) is designed to connect to the F&P Neopuff™ Infant Resuscitator (RD900 Series) or a gas-powered resuscitator compliant to ISO 10651-5:2006.

Product Specifications

Input Flow Range	5-15 L/min
Circuit Length	1.6 m (63")
Positive End-Expiratory (PEEP)*	@5 L/min 1 to 5 cmH ₂ O [mbar] @8 L/min 1 to 9 cmH ₂ O [mbar] @10 L/min 2 to 15 cmH ₂ O [mbar] @15 L/min 3 to 25 cmH ₂ O [mbar]
Device Connection	10 mm Female/15 mm Male
T-Piece Connection	15 mm Female and 22 mm Male (Conical Connectors)
Disposal	Discard according to hospital protocol
Environmental limits	Operating: -18 °C to 50 °C (-0.4 to +122 °F), up to 95% humidity Storage: -10 °C to 50 °C (+14 to +122 °F), up to 95% humidity
Patient Mass Range	<10 kg (22 lb)

*Values are representative only.

Warnings

⚠️ DO NOT connect the T-Piece Circuit (RD1300-10) directly to a high pressure gas source. The T-Piece Circuit (RD1300-10) must be connected to the outlet port of the Neopuff™ or gas-powered resuscitator compliant to ISO 10651-5:2006.

⚠️ Bypassing the resuscitation device may result in serious patient harm (e.g. lung injury).

Ⓜ️ The T-Piece Circuit is designed for SINGLE PATIENT USE ONLY. Reuse may result in transmission of infectious substances, interruption to treatment, or serious harm.

⚠️ DO NOT soak, wash, sterilize or re-use this product.

⚠️ DO NOT occlude the resuscitation tubing.

⚠️ Ensure no smoking, naked flames or sources of ignition are present when using this circuit.

⚠️ The T-Piece circuit is only to be used with a flow rate from 5-15 L/min. DO NOT attempt to use product outside this flow range.

⚠️ The T-Piece circuit and Neopuff™ are intended for use by adequately trained medical practitioners.

⚠️ California residents please be advised of the following, pursuant to Proposition 65:

This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm. For more information, please visit: www.fphcare.com/prop65

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