
F&P WigglewING™

Next generation care


FOR USE WITH
F&P Optiflow™ Junior 2
F&P Optiflow™ Junior 2+

INSTRUCTIONS FOR USE


I F&P WigglewING™ Fitting


A. RIP into two

B. STICK to face



C. FOLD to secure NG tube


REF 607831 REV C 2021-03 © 2021 Fisher & Paykel Healthcare Ltd

Manufacturer ■ Fisher & Paykel Healthcare Ltd, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348 Panmure, Auckland 1741, New Zealand Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.co.nz
Web: www.fphcare.com Importer/ Distributor **Australia** (Sponsor) Fisher & Paykel Healthcare Pty Ltd, 19-31 King Street, Nunawading, Melbourne, Victoria 3131. Tel: +61 3 9871 4900 Fax: +61 3 9871 4998 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22
Benelux Tel: +31 40 216 3555 Fax: +31 40 216 3554 **Brazil** Fisher & Paykel do Brasil, Rua Sampaio Viana, 277 cj 21, Paraíso, 04004-000, São Paulo – SP, Brazil Tel: +55 11 2548 7002 **China** 代理人/售后服务机构: 费雪派克医疗保健 (广州) 有限公司, 广州高新技术产业开发区科学城科丰路31号G12栋301号 电话: +86 20 32053486 传真: +86 20 32052132
Denmark Tel: +45 70 26 37 70 Fax: +46 83 66 310 **Finland** Tel: +358 9 251 66 123 Fax: +46 83 66 310 **France** Fisher & Paykel Healthcare SAS, 10 Av. du Québec, Bât F5, BP 512, Villebon-sur-Yvette, 91946 Courtaboeuf Cedex, France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 Email: c.s@fphcare.fr **Germany** Fisher & Paykel Healthcare GmbH, Deutschland, Österreich, Schweiz, Wiesenstrasse 49, D73614 Schorndorf, Germany Tel: +49 7181 985990 Fax: +49 7181 985996
Hong Kong Tel: +852 2116 0032 Fax: +852 2116 0085 **India** Tel: +91 80 2309 6400 **Ireland** Tel: 1800 409 011 Fax: +44 1628 626 146 **Italy** Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 5117 7110 Fax: +81 3 5117 7115
Korea Tel: +82 2 6205 6900 Fax: +82 2 6309 6901 **Mexico** Tel: +52 55 9130 1626 **Norway** Tel: +47 21 60 13 53 Fax: +47 22 99 60 10 **Poland** Tel: 08000 49 119 77 Fax: 08000 49 119 78 **Russia** Tel. and Fax: +7 495 782 21 50
Spain Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8564 76680 Fax: +46 8366310 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **Taiwan** Tel: +886 2 8751 1739 Fax: +886 2 8751 5625 **Turkey** İthalatçı Firma: Fisher Paykel Sağlık Ürünleri Ticaret Limited Şirketi, İletişim Bilgileri: Östım Mahallesi 1249. Cadde No:6, Yenimahalle, Ankara, Türkiye 06374, Tel: +90 312 354 34 12 Fax: +90 312 354 31 01 **UK** Fisher & Paykel Healthcare Ltd, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, UK Tel: 0800 132 189 Fax: +44 1628 626 146 **USA/Canada** Tel: 1800 446 3908 or +1 949 453 4000 Fax: +1 949 453 4001.

F&P WigglewiNG™ Removal

F&P Optiflow™ Junior 2 ¹ Cannula Size	F&P WigglewiNG
 OJR410 XS	WJR210 
 OJR412 S	
 OJR414 M	WJR212 
 OJR416 L	
 OJR418 XL	

F&P Optiflow™ Junior 2+ ¹ Cannula Size	F&P WigglewiNG
 OJR520 XXL	WJR214 

1. Refer to the F&P Optiflow™ Junior 2 and F&P Optiflow™ Junior 2+ nasal cannula user instructions for appropriate sizing.

F&P WIGGLEWING NASAL CANNULA

English en

INTENDED USE

The F&P WigglewiNG is intended to secure nasogastric tubing to the face in conjunction with the F&P Optiflow Junior 2 and F&P Optiflow Junior 2+ nasal cannula ranges. The F&P WigglewiNG is a single use accessory intended for use on neonates, infants and children in a hospital environment.

CONTRAINDICATIONS/SIDE EFFECTS/CLINICAL BENEFITS

For Contraindications, Side Effects & Clinical Benefit Refer to Interface User Instruction

GENERAL WARNINGS

- This product is only designed and verified for use with equipment, accessories and spare parts approved by F&P. Unauthorized equipment, accessories or spare parts which are used with this product may impair performance of this product or compromise safety (including potentially causing serious patient injury).
- Appropriate patient monitoring (e.g. oxygen saturation) must be used at all times. Failure to monitor the patient (e.g. in the event of an interruption to gas flow) may result in serious injury or death.
- If using supplemental oxygen, keep ignition sources away from the patient.
- Pre-existing craniofacial abnormalities, deformities, malformations or trauma may be exacerbated by NHF interfaces and/or the retention mechanism and may not permit therapy to be delivered as intended, leading to further injury or death.
- The use of this product is not without risk, even if used as intended. Following all instructions and warnings provided, risks of hypoxic injury and skin damage remain. These risks may result in serious injury or death.
- Application of NHF is known to generate positive airway pressure which may exacerbate preexisting, untreated serious air leak syndrome which may lead to further serious injury or death.
- Any serious incident that has occurred in relation to this device should be reported to your Fisher & Paykel Healthcare representative and the local competent authority.

GENERAL CAUTIONS

- This product is intended to be used for a maximum of 7 days.
- Regularly monitor the patient to ensure skin integrity and that the skin underneath the F&P WigglewiNG remains dry. A barrier film may be used between the NG tubing and the patient's upper lip to prevent irritation.
- DO NOT soak, sterilize or reuse this product. Avoid contact with chemicals, cleaning agents, or hand sanitizers.
- Reuse may result in the transmission of infectious substances, interruption to treatment, serious injury or death.
- DO NOT place the F&P WigglewiNG on the patient's eyes, ears or injured skin.
- DO NOT stretch the cannula on application; this may cause increased pressure to the patient's skin. If necessary, the cannula may be repositioned.
- Ensure the cannula is placed directly on the F&P WigglewiNG. Direct skin contact caused by cannula misalignment may result in skin breakdown.
- DO NOT use if the product or its packaging has been tampered with.
- Failure to apply and use this product within the directions, transport, storage and operating conditions specified in the labeling and user instructions may impair performance of this product or compromise safety (including potentially causing serious patient injury).

F&P, Optiflow, Wigglepads, WigglewiNG, and the Seahorse, Crab, Starfish, Octopus, Turtle, Dolphin and Stingray images shown in this user instruction are trade marks of Fisher & Paykel Healthcare Limited.

© 2020 Fisher & Paykel Healthcare Limited.

For patent information, see www.fphcare.com/ip

California residents please be advised of the following, pursuant to Proposition 65:

WARNING

This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm.

For more information, please visit: www.fphcare.com/prop65

I F&P WIGGLEWING FITTING

For the correct F&P WigglewiNG positioning, refer to the placement of the F&P Wigglepads™ 2 in the F&P Optiflow™ Junior 2 and F&P Optiflow Junior 2+ nasal cannula user instructions.

A. RIP INTO TWO

R PATIENT'S RIGHT NOSTRIL L PATIENT'S LEFT NOSTRIL

- Prepare the patient's skin according to hospital protocol.
1. Determine which nostril you would like to place the NG tube in.
 - Hold the F&P WigglewiNG over the patient's cheek on the same side as this nostril.
 2. Face the white tabs (marked **2**) towards you and the hook pads away from you.
 - Place the straighter edge towards the patient's nose.
 3. Rip off the entire upper third of the F&P WigglewiNG (distal to the patient's chin) by holding the middle section and ripping off the upper section.
 4. This will leave you with the lower 2/3 of the F&P WigglewiNG.

B. STICK TO FACE

5. Peel off the blue tab (marked **1**).
6. Stick the F&P WigglewiNG to the patient's cheek
7. Peel off the white tabs (marked **2** and F&P WigglewiNG).
 - Insert the NG tube according to hospital protocol.
8. Place the NG tube along the perforations.

C. FOLD TO SECURE THE NG TUBE

9. **Fold** the F&P WigglewiNG over to secure the NG tube.
10. Apply the section removed in Step 3 as you would a Wigglepad on the opposite cheek.
 - Remove the existing Wigglepad from the nasal cannula.
11. Secure the cannula in place so the cannula bridge rests close to the nose without touching the septum. Do not stretch the cannula during application.

II F&P WIGGLEWING™ REMOVAL

- Place fingertip on the outside edge of the F&P WigglewiNG and gently peel nasal cannula away. Starting from the outside, peel towards the nose.
1. Lift the edge of the F&P WigglewiNG. Use a damp cloth to wipe the patient's skin and the underside of the F&P WigglewiNG while gently peeling away from the patient's face.
 2. Holding the NG tube securely between the patient's nose and the F&P WigglewiNG, detach the NG tube from the F&P WigglewiNG by tearing along the perforations.

SYMBOL DEFINITIONS

	Do not reuse		This product is not made with natural rubber latex
	Not made with phthalates (DEHP, DBP, BBP)		Storage temperature range
	Batch code		Catalogue number
	Consult instructions for use	Rx only	Prescription only
	Authorized representative in the European Community	CE 0123	CE Mark
	Use by date		Medical device
	Date of manufacture		Manufacturer
	Recyclable PET		Caution
	Maximum 7 days use		

DISPOSAL

There are no specific disposal requirements for this device. Disposal protocols appropriate to single use devices that may be contaminated with potentially infectious substances should be followed.

CAUTION

Using this product beyond 7 days may impair performance of this product or compromise safety (including potentially causing serious patient injury).

OPERATING CONDITIONS

Ambient temperature range: 18 to 26 °C.

NOTE: This product may be used in an incubator.

APPROVED F&P RESPIRATORY INTERFACES

F&P Optiflow™ Junior 2 nasal cannula

F&P Optiflow™ Junior 2+ nasal cannula

Refer to the F&P Optiflow Junior 2 and F&P Optiflow Junior 2+ nasal cannula user instructions for appropriate sizing.

California residents please be advised of the following, pursuant to Proposition 65:

WARNING This product contains chemicals known to the State of California to cause cancer, birth defects and other reproductive harm.

For more information, please visit: www.fphcare.com/prop65

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