





<b>WigglewiNG</b> 	<b>WJR210</b>		<b>WJR212</b>			<b>WJR214</b>
						
Weight (per ea)	1.05 g		1.87 g			3.38 g
Interface size	XS	S	M	L	XL	XXL
Description	The F&P WigglewiNG is intended to secure nasogastric (NG) tubing to the face in conjunction with the F&P Optiflow Junior 2 and F&P Optiflow Junior 2+ nasal interfaces.					

**PRODUCT**

Quantity	<b>Case:</b> 20 each (2 packs of 10)	<b>Case:</b> 10 each (1 pack)
Box components	F&P WigglewiNG, user instructions	
Packaging materials	Outer: cardboard box; user instruction: paper	
Packaging dimensions (L x W x H mm)	<b>Pack:</b> 113 x 82 x 25; <b>Case:</b> 171 x 123 x 32	<b>Case:</b> 171 x 123 x 32

**PERFORMANCE**

Ambient range	18 – 26 °C / 64 – 79 °F
Usage	Single patient use
Duration of use	7 days
Storage temperature	<b>Min:</b> -10 °C / 14 °F, <b>Max:</b> +50 °C / 122 °F
Shelf life	3 years

**COMPONENTS AND COMPOSITIONS**

Predominant materials	Thermoplastic elastomer, hydrocolloid
Materials not present	Not manufactured with natural latex, PVC or phthalates (DEHP, DBP, BBP)
Manufacturing mode	Non-invasive device; produced in a controlled working environment
Disposal	Incineration or according to hospital protocol

**REGULATORY**

Classification	AU-IIa, EU-IIa, Canada-IIa, US-I	
Country of origin	New Zealand	
Notified body	TÜV SÜD Product Services GmbH CE0123	
GTIN Codes - WigglewiNG Case	<b>WJR210</b>	09420012454836
	<b>WJR212</b>	09420012454843
	<b>WJR214</b>	09420012454867
GMDN Code	16866	